



COVID-19 Pandemic: Issues with Laboratory Tests when Authorized for Emergency Use (EUA)

Presented by


Paulo Pereira, Ph.D., Senior Researcher, R&D, Portuguese Institute of Blood and Transplantation, Lisbon, Portugal

Paulo Pereira, Ph.D.

- Doctor of Philosophy, Biotechnology, specialization Microbiology
- Senior Researcher, Head of the R&D Department, Portuguese Institute of Blood and Transplantation, Lisbon, Portugal
- Contributing Editor in Westgard QC
- Technical Expert on CLSI Document Development Committee on EP12-A3
- Author of several publications on quality control in the medical laboratory and in the Blood, Cells, Organs and Tissues Bank
- Author of the book "Quality control of qualitative tests for medical laboratories" (2019)



Introduction to the SARS-Cov-2 virus and COVID-19



World Health Organization

Search by Country, Territory, or Area

COVID-19

China / Society

Covid-10 Testes serológicos

Quote of the day:

"[T]here are many tests on the market and many of them - now that I can speak very freely that I'm not in the government - many of them are junk."


Dr. Tom Frieden, former director US CDC. Interviewed on STAT's Helen Branswell.

Hospital de Santa Maria

emir 21.04.2020 às 8h03

cent strike rate

WESTGARD QC



World Health Organization

Search by Country, Territory, or Area

COVID-19

China / Society

Covid-10 Testes serológicos

Quote of the day:

"[T]here are many tests on the market and many of them - now that I can speak very freely that I'm not in the government - many of them are junk."

Dr. Tom Frieden, former director US CDC. Interviewed on STAT's Helen Branswell.

Hospital de Santa Maria


emir 21.04.2020 às 8h03


cent strike rate

WESTGARD QC




World Health Organization







What is happening with standardization and good practices?



What is happening with standardization and good practices?


- Europe and America were not prepared for an outbreak
- We have not learned from the lessons of Asia and Africa
- We have not learned from the history of pandemics
- Regulatory agencies skipped routine IVD assessments for SARS-Cov-2 tests, so they could be available quickly in the medical laboratory
- This approach suppressed strict supervision, which limits, from a clinical perspective, fitness for purpose / clinical decision
- Some manufacturers are still developing tests with the same scientific rigor, others may not be






What is happening with standardization and good practices?

- It is crucial to gradually implement the prequalifications required of IVD manufacturers, including specifications for sampling and variability between reagent lots
- The risk to clinical decisions based on laboratory results is higher
- Clinical decisions must be made incorporating many other variables, contrary to what would be expected with a screening test with high clinical sensitivity
- **Good laboratory practices have been reduced**
- **The risk of false results is increased**







What is happening with standardization and good practices?

• Process approach


| Pre-Pre-Examination Processes | Pre-Examination Processes | Examination Processes | Post-Examination Processes | Post-Post-Examination Processes |
|--|--|--|---|---|
| <ul style="list-style-type: none"> • Good practices • The SARS-Cov-2 RNA RT-PCR should be selected for screening • The serological test for IgG / IgM tests should be selected to assess immunity | <ul style="list-style-type: none"> • Good practices • Non-compliance with causes in training | <ul style="list-style-type: none"> • Good practices • Quality of examination results must be ensured | <ul style="list-style-type: none"> • Good practices • Non-conformities in the results reports | <ul style="list-style-type: none"> • Good practices • Given the risk of false negatives, clinical decisions should also match the symptoms and risk group • Window period should be considered |






What is happening with standardization and good practices?


- Several guidelines published by the WHO to fast-track testing
- Minimal harmonization of good laboratory practices
- It aims to get a quick response from the reference laboratories (“in-house” tests) and the IVD manufacturers






What is happening with standardization and good practices?

- WHO interim guidance for laboratory testing ([updated March 19, 2020](#)) and laboratory testing strategy recommendations for COVID-19 ([updated March 22, 2020](#))
- Molecular assays to diagnose COVID-19 ([FIND Web site](#) and [shared protocols for “in-house” developed molecular assays](#))
- WHO reference laboratories providing confirmatory testing for COVID-19 (WHO reference laboratories providing confirmatory testing for COVID-19, [updated March 2](#), and booking form for national laboratories, [updated March 11, 2020](#))
- Guidance for laboratories shipping specimens to WHO reference laboratories that provide confirmatory testing for COVID-19 virus ([updated March 31, 2020](#))
- WHO interim guidance for laboratory biosafety related to COVID-19 virus ([updated March 22, 2020](#))





What are the pros of “Emergency Use Only” validation?



What are the pros of “Emergency Use Only” validation?


- Pros
 - Simplification of the validation methodology
 - Low cost
 - Fast process of designing and developing a new kit
 - Rapid implementation of new kits on the market
 - Rapid reporting of laboratory results
 - Clinical decisions based on laboratory results








What are the pros of “Emergency Use Only” validation?


- Potential for...
 - Less quality and regulation (poor testing)
 - Recognized manufacturers are developing new tests in a stressed production scenario
 - Manufacturers of unknown quality are producing reagent kits
 - Supply chains are strained, presenting a high risk of shortage at several points, such as the availability of raw materials and reagent kits on the market
 - The effect of false results in an epidemic outbreak has been misestimated






What are the pros of “Emergency Use Only” validation?




What are the pros of “Emergency Use Only” validation?



Matrix Skills Good Practices Standards/ Harmonization Regulatory Compliance Healthcare System

QC


What is going on with *in vitro* diagnostic (IVD) medical devices?



What is going on with *in vitro* diagnostic (IVD) medical devices?

- Fast answer to the needs of tests by the manufacturers
- Introduction of new nucleic acid amplification tests (NAAT)
- Later introduction of serology, at this moment, mainly POCT
- The validation reports are a very simplified version and not harmonized to the commonly required European Commission or the FDA requirements
- Emergence of new manufacturers
- The performance verification in reference laboratories of some tests of new manufacturers led to their rejection
- Business opportunity

QC





What is going on with *in vitro* diagnostic (IVD) medical devices?


- Balance

Reliable and consistent QC practices

Fast availability of new tests








What is going on with *in vitro* diagnostic (IVD) medical devices?


- SARS-Cov-2 molecular assay independent evaluation by WHO/FIND/University Hospitals of Geneva

| Manufacturer | Gene target | LoD* | Clinical sensitivity | Clinical specificity |
|--------------|-------------|--------|------------------------|------------------------|
| A | E | 1-10 | 92% (95% CI: 81, 97) | 100% (95% CI: 96, 100) |
| A | S | 1-10 | 100% (95% CI: 96, 100) | 100% (95% CI: 96, 100) |
| B | ORF1 | 1-10 | 100% (95% CI: 93, 100) | 99% (95% CI: 95, 100) |
| C | E | 10-50 | 100% (95% CI: 93, 100) | 100% (95% CI: 96, 100) |
| C | RdRP | 50-100 | 90% (95% CI: 79, 96) | 98% (95% CI: 93, 99) |
| D | S | 1-10 | 100% (95% CI: 93, 100) | 100% (95% CI: 96, 100) |
| E | RdRP | 10-50 | 100% (95% CI: 93, 100) | 100% (95% CI: 96, 100) |
| F | E | 1-10 | 100% (95% CI: 96, 100) | 100% (95% CI: 96, 100) |

* Copies per reaction



What is the window period for NAAT and serology tests?




What QC suggestions can be given in this outbreak?

What QC suggestions can be given in this outbreak?

- Statistical tools in the outbreak
- Binary examination agreement
- Clinical performance evaluation
- Evaluation of the limit of detection (LoD) in NAAT, i.e., minimum detectable concentration
- Cross-reactivity studies
- Serotypes/Genotype variation
- Stability

QC




What QC suggestions can be given in this outbreak?


- Binary examination agreement

| Candidate test results | Comparative test results | | Total |
|------------------------|--------------------------|----------|-------|
| | Positive | Negative | |
| Positive | a | b | a + b |
| Negative | c | d | c + d |
| Total | a + c | b + d | N |

A free tool is available at the Westgard QC Web site <https://www.westgard.com/> to compute binary examination agreement

- PPA[%] = $a / (a + c) * 100$
- PNA[%] = $d / (b + d) * 100$





What QC suggestions can be given in this outbreak?


- Binary examination agreement


| Candidate test results | Comparative test results | | Total |
|------------------------|--------------------------|----------|-------|
| | Positive | Negative | |
| Positive | 20 | 2 | 22 |
| Negative | 0 | 18 | 18 |
| Total | 20 | 20 | 40 |

Target PPA = 90%, the target of the 95% CI low limit is 75%

Target PNA = 80%, the target of the 95% CI low limit is 70%


- PPA = $(20 / 22) * 100 = 90\%$ (95% CI: 84% - 100%)
- PNA = $(18 / 20) * 100 = 90\%$ (95% CI: 70% - 97%)






What QC suggestions can be given in this outbreak?

- Binary examination agreement
- The **key point** is to have a **test with very high clinical accuracy as a comparator**
- There will always be an uncertainty in the agreement of the results
- If the comparator is a test with worse clinical accuracy than the new test, the results may be misinterpreted
- If the test is a "gold standard," the clinical accuracy is evaluated and not the agreement
- Associated to rare tests, such as some "in-house" tests






What QC suggestions can be given in this outbreak?


- Clinical performance evaluation

| Candidate test results | Clinical accuracy criteria | | Total |
|------------------------|-----------------------------|-----------------------------|---------|
| | Positive | Negative | |
| Positive | True-positive (TP) results | False-positive (FP) results | TP + FP |
| Negative | False-negative (FN) results | True-negative (TN) results | FN + TN |
| Total | TP + FN | FP + TN | N |

- $Se[\%] = TP / (TP + FN) * 100$
- $Sp[\%] = TN / (FP + TN) * 100$

A spreadsheet can be downloaded from MedLabQuality Web site <https://medlabquality.com> to compute the clinical performance evaluation, and LoD using probit regression





What QC suggestions can be given in this outbreak?


- Clinical performance evaluation


| Candidate test results | Clinical accuracy criteria | | Total |
|------------------------|----------------------------|----------|-------|
| | Positive | Negative | |
| Positive | 30 | 0 | 30 |
| Negative | 0 | 30 | 30 |
| Total | 30 | 30 | 60 |

- Sensitivity = $(30 / 30) * 100 = 100\%$ (95% CI: 89% - 100%)
- Specificity = $(30 / 30) * 100 = 100\%$ (95% CI: 89% - 100%)

Target sensitivity = 100%, the target of the 95% CI low limit is 85%

Target specificity = 90%, the target of the 95% CI low limit is 80%






What QC suggestions can be given in this outbreak?

- Clinical performance evaluation

- A number of samples must be taken to allow **statistical robustness** of calculations: infected individuals (D_1), and healthy individuals (D_0)
- The **clinical robustness** of the evaluation requires samples that are epidemiologically representative of the tested population
- The 95% CI should be used, since it illustrates the statistical robustness, as well as allowing the inference for a population with **the same characteristics of the samples**
- The number of samples must suit the target of the evaluation
- For example, for a n of 30 for the determination of clinical sensitivity, the 95% CI can never be tighter than the range of 87% to 100%





What QC suggestions can be given in this outbreak?

- Clinical performance evaluation
- 95% CI is related to the **statistical power** of evaluation
- The **clinical power** of evaluation is related to the representativeness of the infected samples
- Samples of infected individuals are available in the outbreak
- Negative sample could be taken from a serotec or rejected plasma bags





What QC suggestions can be given in this outbreak?

- Is it possible to have robust validations during an outbreak?
- Yes, but...
- Due to several limitations, the first validations must be carried out in reference laboratories according to harmonized practices
- The biggest limitation is the limited number of infected samples and their limited variability at the beginning of the outbreak
- All commercial tests should be **periodically revalidated**, mainly by national agencies, and the reported results should be public
- National agencies should purpose performance goals
- Performance targets should be **reviewed periodically**






What QC suggestions can be given in this outbreak?


- Goal 1**
Improvement of good quality control practices
- Goal 2**
Clinical performance evaluation based on clinical sensitivity and clinical specificity
- Goal 3**
Compliance assessment



What QC suggestions can be given in this outbreak?



What QC suggestions can be given in this outbreak?



Conclusion



Conclusion

- There needs to be a plan for outbreaks like this, including for waves
- The importance of a reliable and consistent med lab results should be reinforced
- Regulatory agencies of IVD manufacturers and medical laboratories should also be prepared for this type of outbreak
- The regulation of IVD medical devices should be strengthened
- Preventive actions should be implemented to avoid even more serious consequences in pandemic outbreaks with more pathogenic agents
- WHO's importance and strength should be reinforced





Links and resources

- [Webinar handouts \(MedLabQuality\)](#)
- [COVID-19 spreadsheet \(MedLabQuality\)](#)
- [Binary examination agreement online calculator \(Westgard QC\)](#)
- [Westgard QC Lesson - Basic Validation of Qualitative Tests](#)
- [European Centre for Disease Prevention and Control of European Union](#)
- [FIND Diagnostics - COVID-19 Diagnostics Research Centre \(Data Base of Assays and Independent Test Evaluations\)](#)
- [The International Federation of Clinical Chemistry and Laboratory Medicine \(IFCC\) - Information Guide on COVID-19](#)





Links and resources

- [US FDA- Coronavirus Disease 2019 \(COVID-19\)](#)
- [US Centers for Disease Control and Preventions - Laboratories](#)
- [The College of American Pathologists](#)
- [The Royal College of the Pathologists of Australasia](#)
- [NRL Australia](#)
- [Zhejiang University School of Medicine \(Handbook of COVID-19 Prevention and Treatment\)](#)
- [WHO - Country & Technical Guidance - Coronavirus disease \(COVID-19\)](#)
- [WHO database of publications on COVID-19](#)





COVID-19 Pandemic: Issues with Laboratory Tests when Authorized for Emergency Use (EUA)

Presented by

Paulo Pereira, Ph.D., Senior Researcher, R&D, Portuguese
Institute of Blood and Transplantation, Lisbon, Portugal
